

UNICORN 7.6.0

1 (2)

User: Default 05/02/2026 09:00:39 +00:00

AKTAPure4

System typeÄKTA Chromatography system
Instrument configurationAKTA pure 25 (1.13.0.63)
System isactive
Instrument ServerSB-AKTA04
Instrument serial number2802373
Connection based onInstrument serial number

Instrument Configuration

Name.....AKTA pure 25
Version.....1.13.0.63

AKTA pure 25 instrument configuration

Strategy Name.....AKTApure25_322

. Predefined Method Configuration

Name.....AKTA pure
Version.....1.7.0.2

Phase Configuration

Name.....AKTA pure
Version.....1.7.0.2

Performance Test Configuration

Name.....AKTA pure
Version.....1.7.0.2

Process Picture Configuration

Name.....AKTA avant pure
Version.....1.2.2.123

Customization Configuration

Name.....AKTA pure
Version.....1.0.4.0

VID Configuration

Name.....AKTA pure 25
Version.....13.0.0.91

Firmware Configuration

Name.....AKTA pure
Version.....4.21.20.0

Components.

Injection valve.....V9-Inj
Sample pump S9
Inlet A.....V9-IA (7-ports)
Inlet B.....V9-IAB part B (2-ports)

UNICORN 7.6.0

2 (2)

User: Default 05/02/2026 09:00:39 +00:00

Sample inlet (V9-IS)

Loop valve (V9-L)

V9-C (5-columns)

Outlet valve.....V9-Os (1-outlet)

UV monitor.....U9-L (fixed)

Conductivity monitor (C9)

Fraction collector.....Fraction collector (F9-C)

System pump A (P9 A)

System pump B (P9 B)

System pressure monitor (R9)

Mixer (M9)